# PCR-RFLP directly from blood using Thermo Scientific Phusion Blood Direct PCR Kit

Pia Kuusisto, Thermo Fisher Scientific

#### **Key words**

- RFLP
- Direct PCR
- Genotyping
- Human Blood
- SNP Alleles
- Restriction
  Digest

#### **Abstract**

This application protocol describes how PCR-RFLP (restriction fragment length polymorphism) assays can be performed from blood samples without DNA purification. In the first step, genomic DNA is amplified directly from whole blood using Thermo Scientific Phusion Blood Direct PCR Kit. The PCR product is subsequently digested with a restriction enzyme that recognizes a specific sequence of interest (e.g. a SNP site) in this DNA fragment. Based on the sample genotype, the digestion produces DNA fragments which are then analyzed by agarose gel electrophoresis. Presented here is a specific example of a PCR-RFLP based SNP (single-nucleotide polymorphism) genotyping assay performed directly from blood without DNA purification before PCR or digestion (see Figure 1A for assay design).

#### **Materials and methods**

- Human whole blood with EDTA as anticoagulant
- Phusion® Blood Direct PCR Kit
- HhaI
- Thermo Scientific Piko 24-well Thermal Cycler
- Thermo Scientific Piko UTW (ultra-thin walled) Piko® PCR Plates
- Primers: (Forward and Reverse)
  - 429 bp fragment of human eukaryotic translation initiation factor 2-alpha kinase 3 gene (EIF2AK3)

F:AGCTCCTATAGTAACCTCTTCTTGA ACTCACTTG 34 nt Tm=68.6°C

R: GCTTTCACGGTCTCGGTCCCACTG 24 nt Tm = 75.5°C

Table 1. Reaction conditions for PCR


Component	20 μL reaction	50 μL reaction	Final conc.			
H <sub>2</sub> 0	Add to 20 µL	Add to 50 µL				
2x Phusion Blood PCR Buffer	10 μL	25 μL	1x			
Primer F (Forward)	x μL	x μL	0.5 μΜ			
Primer R (Reverse)	x μL	x μL	0.5 μΜ			
Phusion Blood DNA Polymerase	0.4 μL	1 μL				
Whole blood*	1 μL	2.5 μL				
Optional components for reaction optimization*						
50 mM MgCl <sub>2</sub>	0.6 μL	1.5 μL				
50 mM EDTA	0.5 - 1.0 μL	1.25 - 2.5 μL				
DMS0	1.0 µL	2.5 μL	5%			

<sup>\*</sup> See Phusion Blood Direct PCR Kit manual for more instructions related to optional components.

Table 2. Cycling protocols

	2-step protocol		3-step protocol		
Cycle step	Temp.	Time	Temp.	Time	Cycles
Lysis of cells	98°C	5 min	98°C	5 min	1
Denaturation Annealing* Extension*	98°C - 72°C	1 s - 15-30 s /kb	98°C x°C 72°C	1 s 5 s 15-30 s /kb	35-40
Final extension	72°C 4°C	1 min hold	72°C 4°C	1 min hold	1

<sup>\*</sup> See Phusion Blood Direct PCR Kit manual for more instructions.






Figure 1. Assay design and results of the SNP genotyping example. A. A 429 bp fragment of human EIF2AK3 gene covering the SNP site of interest was first amplified directly from whole blood of different individuals (5 % blood in 50  $\mu$ L reactions, 2-step cycling protocol, 40 cycles). The unpurified PCR products were subsequently digested as described with Hhal restriction enzyme that recognizes only the G allele in the SNP site of interest. B. The resulting fragments revealing the genotype of each individual were analyzed on agarose gel in 188 bp parallel with the untreated PCR products. M, size marker; A, G and A/G correspond to the SNP alleles of each sample.

#### **Restriction enzyme digestion**

After PCR, the reactions were centrifuged at 1000~x~g for 2 minutes and the supernatants were collected for restriction digestion. The digestions were prepared directly in the supernatants with HhaI (0.4 U/ $\mu$ L). The 10  $\mu$ L reactions were incubated for one hour at 37°C. The resulting fragments were analyzed on agarose gel and compared to untreated PCR products.

## Note 1: Restriction digestion in Phusion Blood PCR Buffer

When amplifying DNA directly from whole blood, the PCR product contains blood and PCR derived components that may interfere with the subsequent digestion. The inhibitor tolerance varies from enzyme to enzyme and therefore further optimization of the reaction conditions may be needed. Several restriction enzymes tolerate inhibitors such as salt and blood constituents quite well and the digestion can be performed directly in a PCR reaction after removing most of the blood components by centrifugation. If the restriction enzyme used is not fully active in Phusion Blood PCR buffer, please follow these instructions below:

- 1. The most convenient option is to dilute the PCR reaction before digestion. Usually diluting the PCR product 1:2 in H<sub>2</sub>O helps to dilute the buffer components and/or inhibitors in the reaction, allowing the restriction enzyme to perform optimally. If the PCR yield is high, it is possible to dilute even more.
- 2. For some restriction enzymes dilution of the PCR reaction is not enough to provide optimal conditions for efficient digestion. In this case, the PCR product needs to be purified before digestion.
- 3. Additionally, further optimization

of the reaction conditions, such as adjustment of the reaction time and restriction enzyme amount, as well as adding restriction enzyme's unique reaction buffer to the digestion, may be needed.

# Note 2: Restriction digestion producing cohesive-ends for cloning

The presence of Phusion Blood DNA Polymerase in digestion reaction may alter the cohesive ends of the DNA fragments that have been cleaved, affecting subsequent ligation. Therefore, for cloning applications, the digestions producing cohesive ends should be performed on purified PCR products.

### **Technical Support**

In North America (US, Canada) and Central/South America Techservice.genomics@ thermofisher.com +1 (800) 235-9880 +1 (303) 604-9499

In Europe, Middle East, and Africa Techservice.emea.genomics@thermofisher.com (+)44 1372 840410

In Other Countries www.thermoscientific.com/ pcrordering

Copyright © 2011 Thermo Fisher Scientific, Inc. All Rights Reserved.

All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Literature Code: 0013511J01U

The Thermo Scientific Direct PCR approach allows for amplification of DNA directly from various starting materials such as blood, mouse ear and tail tissues, plants, and FFPE tissue samples. For more information about the Direct PCR products and protocols, please visit

www.thermoscientific.com/directpcr

