

AP 110 AP 115 AP 125



INSTRUCTION MANUAL



68X452201 Rev. A 07/08


Thank you for selecting an accumet portable (AP) meter. This portable meter is a microprocessor-based instrument with many user-friendly features, all of which are accessible through the keypad. Please read this manual thoroughly as well as your separate electrode manual before operating your instrument.


TABLE OF CONTENTS

1.	THE	INSTRUMENT	. 1
	1.1	DISPLAY	.2
		KEYPAD	
	1.3	CONNECTORS	.4
2.		TURES	
3.	BAT	TERY INSTALLATION	. 5
4.	ELE	CTRODE USE	. 5
5.	pН	OPERATION	
	5.1	Using setup in pH mode to select options	. 6
	5.2	pH Standardization	. 8
	5.3	pH Measurement	. 8
6.	тV	OPERATION	
	6.1	Using setup in mV mode to select options	. 9
	6.2		
7.	R.m	V OPERATION	
	7.1	Using setup in R.mV mode to select options	10
	7.2	R.mV Standardization	11
	7.3	R.mV Measurement	
8.	ION	OPERATION (AP125 only)	12
	8.1	Using setup in the Ion mode	12
		Ion Standardization	
		Measurement in Ion mode	
9.	MEA	SURING TEMPERATURE	15
	9.1	Temperature Standardization (ATC)	15
	9.2	Temperature Standardization (MTC)	
10		TA STORAGE	
	10.1	Storing Value into Memory	17
	10.2	Recalling Value from Memory	17
		ECIFICATIONS	
		ROR MESSAGES	
13	. RE	PLACEMENT PARTS & ACCESSORIES	20

1. THE INSTRUMENT

Before operating the meter, please familiarize yourself with the location and function of its various display elements, keypad and connectors.





1.1 DISPLAY

Figure 2 illustrates the liquid crystal display area. Each of the display elements is referenced in the figure and described below.




Figure 2 Display annunciators

SYMBOL	DESCRIPTION
	Battery indicator with scale
-	Meter operating with optional AC power supply
Setup	Set Up mode
View	View data stored in memory
Std	In Standardization mode
Meas	In Measurement mode
рН	pH Measurement mode
R.	Relative millivolt Measurement mode
mV	Absolute millivolt Measurement mode
ION	Ion Measurement mode
STABLE	Reading stability indicator
EFF%	Slope efficiency in percentage
Clear	Clear (erase) data
AUTO	Auto-read function ON
MM/DD/YY	Month/Date/Year
°C	Temperature in Degrees Celsius
ATC	Automatic Temperature Compensation
mV	Offset in millivolts

1.2. KEYPAD

Figure 3 illustrates the keypad area. The function of each key is described below.




Figure 3 Meter Keypad

KEY	DESCRIPTION
on/off/ ☆	Press to turn the meter ON. The meter will default to the last measurement mode used. When ON, press this key to turn the backlight ON or OFF (to conserve batteries, the backlight is automatically turned off after 30 seconds). When ON, press and hold this key for 3 seconds to turn the meter OFF.
mode	Press this key to toggle between pH, mV, Relative mV (R.mV) and Ion (AP125 only) measurement modes.
setup/meas	Press this key to enter Setup menu. Press again to scroll through the Setup options. Measure key also serves to refresh value when Auto-Hold is active.
std/enter	Press this key to activate and confirm the standardization. Also to confirm selection or change being made in Setup mode.
save/▲	Press this key to store measured data into memory (up to 200 data sets). Also to increase value or make selection in the Setup mode e.g. for date and time setting scroll up selection.
view/▼	Press this key to recall and select memory location of stored data. Also to decrease value or make selection in the Setup mode e.g. for date and time setting scroll down selection.

1.3 CONNECTORS

Figure 4 illustrates the Top View showing connectors of the AP meter.

Note that the meter is waterproof only when the blue rubber plugs and/or the appropriate ATC probe are used. If the meter is totally submersed, water may enter the BNC connector. If necessary, dry the connector to avoid corrosion. Also, If an AC adapter is connected, the meter is <u>not waterproof</u>.




Figure 4 View of Connectors

The **AC** jack is used to connect an optional power adapter. If an AC adapter is not used, Seal the jack with the blue rubber plug to maintain a waterproof state. An AC adapter prevents battery drain by removing the battery from the circuit.

For pH, mV, R.mV, and Ion measurement, attach combination electrodes with BNC connectors directly to the middle **BNC jack** on the meter. Non-combination (half cell) electrodes can be used with appropriate adapters sold separately.

For best results, automatic temperature compensation (ATC) is recommended. Insert the accumet waterproof ATC connector into the **ATC jack**. Insert the attached blue rubber plug if an ATC probe is not used to maintain a waterproof state.

2. FEATURES

FEATURE	DESCRIPTION	
Display Resolution	Operator-selectable using setup/meas key, choose 8888.8 for 0.1, 888.88 for 0.01, 88.888 for 0.001	
Auto Read	When enabled, display locks on current value when stable indicator appears. Press meas key to unlock and resume reading	
Battery Saver Backlight	When enabled, meter automatically turns off after 20 minutes of inactivity To conserve battery life, the backlight automatically turns OFF after 30 seconds	
Slope Display	View electrode slope and mV offset for current reading with setup/meas key	
Auto-buffer Recognition	Meter selects from 5 (USA or NIST) or 6 (EURO) pre-programmed standards with temperature correction of buffer value	
Standardi- zation	User selectable 1-5 cal points for USA/NIST or 1-6 cal points for EURO standardization.	

3. BATTERY INSTALLATION

To install the battery, perform the following steps:

- 1. Remove the battery cover by loosening the 3 screws that hold the battery cover at the back of the meter using a Phillips screwdriver.
- 2. Disconnect the old 9V battery and install new one into meter battery compartment.
- Make certain the battery wires are positioned so that they do not interfere with the closing of the battery cover. Battery cover edges may damage the wires.
- 4. Replace the battery cover and tighten the 3 screws to secure the cover in place.

4. ELECTRODE USE

This meter can be directly fitted with combination electrodes, with or without built-in ATC. See **Accessories** section for pH/ATC electrodes as well as adapters for use with half cells.

- Condition your electrode by soaking as needed in electrode storage solution, pH 4 or 7 buffer, or KCl solution. Complete hydration can take anywhere from minutes to overnight for very dry electrodes.
- Twist the metal shorting cap to remove it from the BNC jack of the meter. Install your electrode by twisting it to lock it in place. The shorting cap is designed to hang freely when the BNC jack is occupied.
- Rinse the electrode using clean water and gently shake dry between measurements. Store in electrode storage solution, pH 4, pH 7, or KCl solution. If refillable, keep the fill hole open during measurement but close for longer periods of nonuse and storage. Refill when the level of fill solution recedes below the manufacturer's recommended level.

Refer to the manufacturer's electrode instruction manual for specific details on electrode use.

5. pH OPERATION

5.1 Using **setup** in pH mode to select options

From the pH Measure screen:

The **setup** button is a scroll button which allows you to view and change several operating parameters. While in **setup** you may:

- Press mode to return to the Measure screen without making a change or selection
- Press setup and scroll through the operating parameters in the meter
- Press enter to accept the parameter as displayed or to accept a change made to that parameter and return to the measurement mode.
- 1. To activate or de-activate the AUTO hold read function press setup/meas once. Select AUTO

HOLd 'Yes' or 'No' by using the **save/**▲ or **view/**▼ keys. Press **std/enter** to accept.

- To display the % slope efficiency 'EFF%' and mV offset 'mV' press setup/meas twice.
- To erase an existing pH standardization/ calibration press setup/meas until 'Clear buf std' (Clear buffer standard) appears on the screen. Press std/enter to accept or mode to cancel.
- To change the buffer standard group, press setup/meas until 'bUF USA' appears. Press std/enter accept the desired group (USA, NIST, or EURO) when it is displayed.
- To change the pH resolution press setup/meas until '88888' appears. Press std/enter when the desired resolution is displayed below; 88.888 to select 0.001 (AP115 & AP125 only) 888.88 to select 0.01 8888.8 to select 0.1
- 6. To activate or deactivate the Automatic Shut off function press setup/meas until 'A.Off' appears. Select 'Yes' or 'No' by using the save/▲ or view/▼ keys and press std/enter to confirm your selection. If 'Yes' is selected, the meter automatically turns off after 20 minutes of inactivity. Note: Selecting 'No' can result in total battery exhaustion—once on, the meter remains on until the meter is manually turned off.
- To set meter time and date press setup/meas until 'MM/DD/YY' appears. Use the save/▲ or view/▼ keys to select the appropriate digits. Press std/enter to confirm selection and to move from minutes to hours; or from month to day to year.
- 8. To clear data memory press **setup/meas** until 'Clear dAtA' appears. Select 'Yes' or 'No' by using the **save/**▲ or **view/**▼ keys. Press **std/enter** to confirm selection.

5.2 pH Standardization

For best results, standardizing (calibrating) your meter/ electrode system against certified accurate solutions is strongly recommended due to electrode variations and changes in electrode response over time. Daily standardization with fresh buffers is common.

- 1. Press **mode** until the display indicates the pH mode.
- Immerse the rinsed electrode(s) into a buffer from the chosen group (USA, NIST or EURO). Stir moderately if possible. Buffers can be used in any order, however neutral (pH 7) is ideally first.
- 3. Press **std/enter** to access the Standardize screen. 'Std' is now displayed indicating this mode.
- 4. Allow the pH value to reach a constant reading (little variation over 30 seconds as well as the appearance of the STABLE annunciator). Press std/enter again to complete standardization of the value. The meter automatically recognizes the buffer and returns to the Measure mode. Successfully standardized buffer values remain on the display until cleared—even after power off.
- Repeat steps (2) to (4) with subsequent buffers. The meter will briefly display 'EFF%' (efficiency as % slope) before returning to the Measure mode.

If this value is not between 70-120%, '**SLOPE Error**' (Electrode error) is displayed. The standardization value is not accepted and the instrument returns to the Measure mode.

5.3 pH Measurement

- Press mode until the display indicates the pH mode. Immerse the electrode(s) into the sample solution. Stir moderately if possible.
- 2. Allow the pH value to reach a constant reading (little variation over 30 seconds as well as the

appearance of the **STABLE** annunciator). The pH reading may be recorded at this time.

- a. If **AUTO** is not displayed, the auto read hold function is not active, and the meter will continuously monitor the live pH value.
- b. If AUTO is displayed, the meter will lock the measured value on the screen when the STABLE annunciator appears. To release the locked value and obtain a new reading, press setup/meas. AUTO will flash until a stable reading is obtained and the new value is locked. Note: enabling the AUTO read hold function can lead to erroneous values if the electrode is slow to respond such as testing of viscous or dirty samples, changing temperatures, etc.

6. mV OPERATION

AP meters display absolute mV values—referenced to zero mV. Although, calibration doesn't apply to absolute measurement, a check of the zero value can be performed by connecting the BNC shorting cap and observing the mV reading. See **Relative R.mV** section for offset standardization.

- 6.1 Using setup in mV mode to select options
- 1. Press mode to select mV measurement mode.
- To change the mV resolution press setup/meas until '88888' appears. Press std/enter when the desired resolution is displayed below; 8888.8 to select 0.1 8888 to select 1
- To activate or deactivate the Automatic Shut off function press setup/meas until 'A.Off' appears. Select 'Yes' or 'No' by using the save/▲ or view/▼ keys and press std/enter to confirm your selection. If 'Yes' is selected, the meter automatically turns off after 20 minutes of inactivity. Note: Selecting 'No' can result in total battery

exhaustion—once on, the meter remains on until the meter is manually turned off.

- To set meter time and date press setup/meas until 'MM/DD/YY' appears. Use the save/▲ or view/▼ keys to select the appropriate digits. Press std/enter to confirm selection and to move from minutes to hours; or from month to day to year.
- 5. To clear data memory press **setup/meas** until 'Clear dAtA' appears. Select 'Yes' or 'No' by using the **save/**▲ or **view/**▼ keys. Press **std/enter** to confirm selection.
- 6.2 mV Measurement
- 1. Press mode until the meter displays the mV mode.
- 2. Immerse the electrode in sample solution.
- When the meter senses that the displayed value is stable, STABLE will appear under the displayed value. Record the reading at this time.

Note: The AUTO read function is not applicable in the mV and R.mV measurement modes.

7. R.mV OPERATION

AP meters can display relative mV (R.mV) values i.e. all values are referenced to a specified offset mV.

- 7.1 Using setup in R.mV mode to select options
- 1. Press mode to select R.mV measurement mode.
- To view saved offset mV value, press setup/meas once. If no offset exists, '----' will be displayed. Press std/enter to accept the offset value.
- To erase an existing Relative mV standard press setup/meas until 'Clear rEL Std' (Clear relative mV standard) appears. Press std/enter to accept or mode to cancel.

- To change the R.mV resolution press setup/meas until '88888' appears. Press std/enter when the desired resolution is displayed below; 8888.8 to select 0.1 8888 to select 1
- 5. To activate or deactivate the Automatic Shut off function press setup/meas until 'A.Off' appears. Select 'Yes' or 'No' by using the save/▲ or view/▼ keys and press std/enter to confirm your selection. If 'Yes' is selected, the meter automatically turns off after 20 minutes of inactivity. Note: Selecting 'No' can result in total battery exhaustion—once on, the meter remains on until the meter is manually turned off.
- To set meter time and date press setup/meas until 'MM/DD/YY' appears. Use the save/▲ or view/▼ keys to select the appropriate digits. Press std/enter to confirm selection and to move from minutes to hours; or from month to day to year.
- 5. To clear data memory press **setup/meas** until 'Clear dAtA' appears. Select 'Yes' or 'No' by using the **save/**▲ or **view/**▼ keys. Press **std/enter** to confirm selection.

7.2 R.mV Standardization

- 1. Press **mode** until the display indicates the R.mV mode.
- Press std/enter to access the Standardize screen. Press std/enter again to confirm the displayed offset R.mV value to zero.
- 7.3 R.mV Measurement
- 1. Press mode until the meter displays the mV mode.
- 2. Immerse the electrode in sample solution.

 When the meter senses that the displayed value is stable, STABLE will appear under the displayed value. Record the reading at this time. This value is the millivolt reading relative to the applied offset.

Note: The AUTO read function is not applicable in the mV and R.mV measurement modes.

8. ION OPERATION (AP125 only)

8.1 Using setup in the Ion mode

From the Measure screen:

The **setup** button is a scroll button which allows you to view and change several operating parameters. While in **setup** you may:

- Press **mode** to return to the Measure screen without making a change or selection.
- Press setup and scroll through the operating parameters in the meter.
- Press enter to accept the parameter as displayed or to accept a change made to that parameter.
- To activate or de-activate the AUTO hold read function press setup/meas once. Select AUTO HOLd 'Yes' or 'No' by using the save/▲ or view/▼ keys. Press std/enter to accept.
- 2. To display the slope value in mV/decade, press setup/meas twice.
- To erase an existing lon standardization/ calibration press setup/meas until 'Clear ION Std' (Clear ion standard) appears on the screen. Press std/enter to accept or mode to cancel.
- To change the lon resolution press setup/meas until '888' appears. Press std/enter when the desired resolution is displayed below; 888 to select 123 88 to select 120 8 to select 100

- 5. To activate or deactivate the Automatic Shut off function press setup/meas until 'A.Off' appears. Select 'Yes' or 'No' by using the save/▲ or view/▼ keys and press std/enter to confirm your selection. If 'Yes' is selected, the meter automatically turns off after 20 minutes of inactivity. Note: Selecting 'No' can result in total battery exhaustion—once on, the meter remains on until the meter is manually turned off.
- To set meter time and date press setup/meas until 'MM/DD/YY' appears. Use the save/▲ or view/▼ keys to select the appropriate digits. Press std/enter to confirm selection and to move from minutes to hours; or from month to day to year.
- 7. To clear data memory press **setup/meas** until 'Clear dAtA' appears. Select 'Yes' or 'No' by using the **save/**▲ or **view/**▼ keys. Press **std/enter** to confirm selection.

8.2 Ion Standardization

For best results, standardizing (calibrating) your meter/ electrode system against certified accurate solutions is strongly recommended due to electrode variations and changes in electrode response over time. Daily standardization with fresh buffers is common.

Please refer to your ion selective electrode manual for calibration standard preparation, maintenance, lonic Strength Adjuster (ISA), storage and other details not provided here.

- 1. Press **mode** until the display indicates the lon mode.
- 2. Press std/enter to access the Standardize screen.
- 3. Press setup/meas to select the desired standard: 0.1, 0.5, 1, 2, 5, 10, 50, 100, 500,1000

- Immerse the ion selective electrode(s) into the standard solution selected in Step (3) adding ISA as needed. Provide moderate stirring if possible.
- 5. Allow the mV value to reach a constant reading (little variation over 30 seconds as well as the appearance of the STABLE annunciator). Press std/enter again to complete standardization of the value. The meter automatically recognizes the standard and returns to the Measure mode. Successfully standardized buffer values remain on the display until cleared—even after power off.
- Repeat steps (2) to (5) with subsequent standards. The meter will briefly display slope value (in mV/decade) before returning to the Measure mode. If this value is not between 15 – 90 mV/decade, 'SLOPE Error' (Electrode error) is displayed. The standardization value is not accepted and the instrument returns to the Measure mode.
- 8.3 Measurement in Ion mode
- Press mode until the display indicates the lon mode. Immerse the electrode(s) into the sample solution. Stir moderately if possible.
- Allow the lon value to reach a constant reading (little variation over 30 seconds as well as the appearance of the STABLE annunciator). The lon reading may be recorded at this time. Note: '----' indicates a two point standardization

Note: '----' indicates a two point standardization has not been completed and the lon reading can not be determined. See lon Standardization.

If **AUTO** is not displayed, the auto read hold function is not active, and the meter will continuously monitor the live pH value.

If **AUTO** is displayed, the meter will lock the measured value on the screen when the STABLE annunciator appears. To release the locked value and obtain a new reading, press **setup/meas**. **AUTO** will flash until a stable reading is obtained and the new value is locked.

Note: enabling the AUTO read hold function can lead to erroneous values if the electrode is slow to respond such as testing of viscous or dirty samples, changing temperatures, etc.

9. MEASURING TEMPERATURE

AP meters can utilize automatic temperature compensation (ATC) or manual temperature compensation (MTC). The ATC may be separate or integrated into a pH electrode. When ATC is used, temperature is continuously updated and visible on the display. When ATC is not used the meter will display 25°C as the default temperature. AP meters can retain simultaneous temperature adjustments for both ATC (with probe) and MTC (default temp when probe is not used). Temperature is displayed in each Measurement mode—pH, mV, R.mV, and Ion.

9.1 Temperature Standardization (ATC)

The temperature reading of your ATC probe can be adjusted at a single temperature value to ensure optimal accuracy. Standardization of ATC is only recommended when temperature errors are suspected and/or when a replacement temperature probe is used.

- 1. Press mode to select mV measurement mode.
- 2. Press **std/enter** to access the Standardize screen. 'Std' is now displayed indicating this mode.
- Dip the ATC probe (or pH/ATC electrode) into a solution of known temperature (i.e. a certified temperature bath). Allow adequate time for the temperature probe reading to stabilize.
- Set the temperature by using the save/▲ or view/▼ keys. The meter allows adjustment up to ±5.0 degrees compared to the factory default.

9.2 Temperature Standardization (MTC)

When not utilizing an ATC probe for best results, the default temperature reading can be adjusted from the factory default value of 25 °C to ensure optimal accuracy. For example, if your pH buffers and samples are at 20 °C and you do not have an ATC probe, it would be desirable to change the MTC from 25 °C to 20 °C. If **'ATC'** is not displayed, MTC is active.

- 1. Press mode to select mV measurement mode.
- 2. Press **std/enter** to access the Standardize screen. 'Std' is now displayed indicating this mode.
- Set the temperature by using the save/▲ or view/▼ keys. The meter allows adjustment to any value between 0.0 to 100.0 degrees.

10. DATA STORAGE

10.1 Storing Value into Memory

The AP 110/115/125 meter can store up to 200 data sets in its non-volatile memory.

In any Measurement mode, press **save/** to store the measured value. "Lo XX Stored (Location number stored) appears on the screen. The value is now stored in the meter's non-volatile memory and the meter returns to Measurement mode.

If memory is full, the meter will prompt whether to overwrite the data (from the first memory location) or not.

To overwrite the data, select 'Yes', then press **std/enter**. The value is now stored in the meter's non-volatile memory and the meter returns to Measurement mode.

If you do not wish to overwrite the data and wish to escape from this menu, select 'No', then press **std/enter**. The meter will return to Measurement mode without storing the data.

10.2 Recalling Value from Memory

In any Measurement mode, press view/▼. The screen will display the latest stored data.

Data retrieval is based on a 'last-in-first-out' basis. To view the specific data; press either **save/** \blacktriangle or **view/** \blacktriangledown keys to select the appropriate data location number. Press **std/enter** to toggle between data and date & time. Press **mode** to return to the Measurement mode.

11. SPECIFICATIONS

рH	AP 110	AP 115	AP 125
Range	-2.00 to		00 to
	20.00	2.00 10	
Resolution	0.1/0.01	0.1 / 0.01 / 0.001	
Relative accuracy	± 0.01	± 0.002	
Input	BNC		
Input Impedance	10 ¹² ohms		
No. of calibration	1 to 6 points		
points	(depending on Buffer selection)		
	USA: 2.00, 4.00, 7.00, 10.00, 12.00		
Buffer values	NIST: 1.	68, 4.01, 6.86,	9.18, 12.45
	EURO: 1.00, 3.00, 6.00, 8.00, 10.00, 13.00		
Min. & max. slope efficiency during calibration	70 – 120%		
mV	AP 110	AP 115	AP 125
mV Range	± 2000		
Resolution	0.1 / 1		
Relative accuracy	±0.2/2		
Relative mV	AP 110	AP 115	AP 125
Absolute mV Range	± 2000		
Resolution	0.1 / 1		
Relative accuracy	±0.2/2		
lon	AP 110	AP 115	AP 125
Range in ppm	''	NO	1x10 ⁻³ to 9.99x10 ⁴
Resolution		NO	1 or 3 digit
Relative Accuracy			0.5% FS
	NO		(monovalent)
			1% FS (divalent)
Input connector		NO	BNC socket
Input Impedance		NO	1 Gohm
No of calibration			2- 5 points
points	1	NO	(min. 2 pts)
Buffer values	1	NO	0.1, 0.5, 1, 2, 5, 10, 50, 100, 500, 1000
Min. & max. slope during calibration	1	NO	15 to 90 mV /decade

Temperature	AP 110	AP 115	AP 125	
Temp. Range °C	-5 to 100			
Temp. Resolution	0.1			
Temp. Accuracy	± 0.3			
Features	AP 110	AP 115	AP 125	
Date & Time	NO	YE	ES	
Auto-Buffer Recognition	YES			
Auto Hold Mode	YES			
Auto Shut Off	Selectable after 20 minutes, Default (on)			
Memory	200 data sets			
Slope/Offset Display	YES			
Backlight	YES			
Ingress Protection	IP 67			
Operating Temperature	0 to 50 °C			
Battery level indicator	YES			
Battery	9 VDC, PP3			
Battery Life	>200 hours without backlight >20 hours with backlight			
Power Adapter	12 VDC			
Dimensions	1.2" x 2.9" x 6.8"		23	
Weight	10 ounces			

12. ERROR MESSAGES

ERROR	DESCRIPTION
	Low battery indicator "LobAt" will show on the LCD before the meter goes off. This indicates that battery power is low and that the batteries need to be replaced.
Or	"Over-range" when reading >20 pH, >2000 mV, or >99900 lon
Ur	"Under-range" when reading <-2 pH or <-2000 mV
CAL Error	Calibration error appears if reading is not within the allowable pH buffer calibration range.
SLOPE Error	Appears if the electrode slope is not within the acceptable range.

13. REPLACEMENT PARTS & ACCESSORIES

Item Description AP110 Meter Only AP110 Meter Kit AP115 Meter Only AP115 Meter Only AP125 Meter Kit AP125 Meter Kit Refillable, pH/ATC electrode, Single Junction Refillable, pH/ATC electrode, Double Junction Gel-filled, pH/ATC electrode, Double Junction Refillable, ORP/Redox electrode Ammonia Ion Selective electrode Chloride Ion Selective electrode Fluoride Ion Selective electrode Sodium Ion Selective electrode	Catalog No. 13-636-AP110A 13-636-AP110 13-636-AP115A 13-636-AP125A 13-636-AP125A 13-636-AP125 13-620-AP50A 13-620-AP50 13-620-AP52 13-620-AP52 13-620-509 13-620-627 13-620-503A
Chloride Ion Selective electrode	13-620-627
Sodium Ion Selective electrode Stainless Steel ATC probe	13-620-629 13-620-503A 13-620-AP53
Hard Portable Meter Carrying Case Kit 110/220 VAC Power Adapter pH 4,7,10 Buffer pack, 500 mL bottles	13-636-AP69 13-636-100 SB105
pH 4,7,10 & Rinse, 20 mL pouches x 5 each	13-300-147

For a complete selection of electrodes and accessories, please refer to the Fisher Scientific Catalog, website, or contact your Fisher Scientific Sales Representative.

To place an order, call 1-800/766-7000, fax 1-800/926-1166, or online www.fishersci.com

For technical support, call 1-888-358-4706 or email accumet@thermofisher.com