

# Get Your pH on Point Calibrate Effectively



**METTLER TOLEDO**

## Preparation



- Switch on the **pH meter** and select the correct buffer group.
- Enter the **correct temperature** for the buffers if no automatic temperature capture is done.
- Pour a sufficient amount of each **buffer solution into a clean beaker**.
- Rinse the electrode** with distilled or deionized water. For refillable electrodes: Ensure that the electrolyte filling hole is open.

## Calibration



- Immerse the electrode into the **first buffer solution** and start the calibration on the pH meter. Start with the lowest pH value.
- Wait until the measurement has **reached the endpoint**.
- Take the electrode out of the buffer solution and **rinse** it.
- Add calibration points by repeating steps a-c with the next buffer solution. Once complete, **end the calibration** on the pH meter.

## Evaluation



- Review the calibration results** on the meter:

	Offset	$\pm (0-20)$ mV	$\pm (20-35)$ mV	$> \pm 35$ mV
<b>Slope</b>				
95–105%				
90–95%				
85–90%				
< 85% or > 105%				

- Electrode in good condition
- Electrode requires cleaning soon
- Electrode requires cleaning and/or regeneration
- Electrode is worn out and needs to be replaced

- Save the calibration** if it is acceptable. The electrode is now **ready** for your measurements.
- If no acceptable calibration can be achieved, visit our **TroubleShooter** and sensor handling movies on: [www.electrodes.net](http://www.electrodes.net)

Visit [eu.fishersci.com](http://eu.fishersci.com) for more information

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2025 Thermo Fisher Scientific Inc. All rights reserved.  
 Trademarks used are owned as indicated at fishersci.com/trademarks.