



For titrating small to medium series where a glass burette would be too slow and a motorised burette too expensive.

- Large LCD clearly shows titrated volume and eliminates meniscus reading errors
- The flowpath of the Digital Burette III is made of highly resistant material: borosilicate glass, PTFE, FEP, PFA, ETFE, platinum-iridium  $Al_2O_3$ -ceramic and PP (screw cap)
- Easy calibration technology, ISO 9000 and GLP compliant
- Long life battery with power for more than 60,000 titrations without recharging
- Enlarged non-slip hand wheels allow titration rate from continuous to drop-by-drop
- Safety discharge system reduces the risk of inadvertent dispensing and splashes if discharge tube is improperly installed or missing
- Available with SafetyPrime™ recirculation valve which avoids reagent waste by recirculating during priming
- Telescopic filling tube adjusts easily to different size bottles, without measuring or cutting
- Zero reset allows multiple titrations
- Accuracy 0.2%, coefficient of variation 0.1%, subdivision 0.01mL
- Supplied with serial numbered performance certificate, conformity certified, operating manual, fine tip titrating tube, telescopic filling tube, valve mounting tool and three polypropylene adapters: GL45/32, GL45/S40 (buttress thread) and GL32/NS 29/32

For accessories and spare parts, please visit our website.

| Catalogue No | Alt. No | Venting valve | Capacity, mL |
|--------------|---------|---------------|--------------|
| BWK-210-011V | 4750150 | Standard      | 25           |
| BWK-210-015N | 4750151 | SafetyPrime™  | 25           |
| BWK-210-031P | 4750160 | Standard      | 50           |
| BWK-210-035H | 4750161 | SafetyPrime™  | 50           |

## Volumetric Solutions

Ready to use solutions for volumetric analysis with a factor up to 0.999. Concentrated solutions are available through the Solutrate range (supplied in singles or as packs of six sealed ampoules).



See the Fisher Scientific Chemicals catalogue for the full range of products.



- Digital burette designed to combine speed and ease of use with high accuracy
- Display rotates through 360° so that it can always be positioned facing the user with the dispensing tube safely pointing away
- Large clear LCD display
- Accuracy  $\pm 0.2\%$ , reproducibility  $\pm 0.01\%$ , subdivision 0.01mL
- Supplied with calibration certificate
- Zero reset for multiple titrations
- Long life battery with power for 60,000 titrations
- Automatic shut off without losing reading to conserve power
- Simple disassembly for easy cleaning with autoclavable liquid pathway
- With BioCote® antimicrobial protection
- 30mm thread supplied with adapters for bottles with 33mm, 38mm and 45mm necks
- Supplied complete with works certificate but without reservoir

| Catalogue No | Alt. No | Capacity, mL |
|--------------|---------|--------------|
| BWK-300-022P | BT30    | 30           |
| BWK-300-042J | BT50    | 50           |

Reservoir base with side neck, amber borosilicate

| Catalogue No | Alt. No | Capacity, mL |
|--------------|---------|--------------|
| BWK-310-020M | DB1400A | 1,400        |
| BWK-310-040G | DB2500A | 2,500        |

Adapters

| Catalogue No | Alt. No | Description     |
|--------------|---------|-----------------|
| BWK-310-501T | BTA40   | For 40mm thread |
| BWK-310-505L | BTA33   | For 33mm thread |
| BWK-310-507H | BTA38   | For 38mm thread |
| BWK-310-509D | BTA45   | For 45mm thread |