

Needles

16



- Stainless steel needles
- Other sizes and point styles available on request

Replacement removable needles (RN) point style 2, bevelled

Catalogue No	Alt. No	Volume, µL	Needle gauge	Needle length, mm	Pack qty
SZR-617-021F	7758-02	0.25 to 100	26s	51	6
SZR-617-031C	7758-03	2.5 to 100	22s	51	6
SZR-617-061Q	7779-02	250 to 10mL	26s	51	6
SZR-617-041W	7779-03	250 to 10mL	22s	51	6
SZR-617-051T	7779-01	250 to 10mL	22	51	6

Point style 3, flat end

Catalogue No	Alt. No	Volume, µL	Needle gauge	Needle length, mm	Pack qty
SZR-617-181G	7768-01	2.5 to 100	26s	51	6
SZR-617-121B	7770-01	2.5 to 100	22s	51	6
SZR-617-191D	7780-01	250 to 10mL	26s	51	6
SZR-617-131V	7780-03	250 to 10mL	22s	51	6
SZR-617-141S	7780-04	250 to 10mL	22	51	6

Kel-F® hub needles for use with LT and TLL syringes

Point style 2, bevelled

Catalogue No	Alt. No	Needle gauge	Needle length, mm	Pack qty
SZR-602-512A	90138	22s	51	6
SZR-602-514T	90126	26	51	6
SZR-602-522U	90122	22	51	6

Point style 3, flat end

Catalogue No	Alt. No	Needle gauge	Needle length, mm	Pack qty
SZR-602-516P	90533	26	51	6
SZR-602-518L	90534	22s	51	6
VAR-930-050L	90134	22	51	6



Syringes, constant rate, 20µL to 200µL

13



For use with gases and liquids.

- Ideal for carbon analysis, waste water analysis, TOC analysis
- Incremental volumes are selected with micrometer precision
- Spring driven plunger injects samples at a constant rate. Field repairable with replacement parts
- Accuracy ±1%
- Sample ejection by push button
- Point style 3, flat end, 51mm needle

Catalogue No	Alt. No	Volume, µL	Increments, µL	Needle gauge
SZR-667-010K	84301	1 to 20	0.05	22s
SZR-667-012R	84303	2 to 50	0.1	22
SZR-667-014N	84302	10 to 200	0.5	22



Reproducibility adapter, Chaney

13



For consistent, reproducible injections.

Adjustable stop prevents plunger damage. Easily installed. For 700, 1700 and 1000 series.

Catalogue No	Alt. No	Volume, µL
SZR-617-500Q	14700	5 to 10
SZR-617-502M	14725	25 to 500

