

HPLC columns, Thermo Scientific Hypersil™ BetaBasic

16



NEW

C₈ columns

Catalogue No	Alt. No	I.D., mm	Particle size, µm	Length, mm
HSA-731-030G	71403-052130	2.1	3	50
HSA-731-035T	71403-152130	2.1	3	150
HSA-731-020J	71403-053030	3	3	50
HSA-731-025W	71403-153030	3	3	150
HSA-731-010M	71403-054030	4	3	50
HSA-731-015C	71403-154030	4	3	150
HSA-730-300J	71403-054630	4.6	3	50
HSA-731-005F	71403-154630	4.6	3	150
HSA-731-085E	71405-052130	2.1	5	50
HSA-731-090L	71405-102130	2.1	5	100
HSA-731-075H	71405-103030	3	5	100
HSA-731-080X	71405-153030	3	5	150
HSA-731-060U	71405-104030	4	5	100
HSA-731-065K	71405-154030	4	5	150
HSA-731-045Q	71405-104630	4.6	5	100
HSA-731-050A	71405-154630	4.6	5	150

C₈ guard cartridges

Catalogue No	Alt. No	I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-731-610L	71403-012101	2.1	3	10	4
HSA-731-605E	71403-013001	3	3	10	4
HSA-731-600X	71403-014001	4 and 4.6	3	10	4
HSA-731-625V	71405-012101	2.1	5	10	4
HSA-731-620Y	71405-013001	3	5	10	4
HSA-731-615B	71405-014001	4 and 4.6	5	10	4

HPLC Columns, Thermo Scientific HyPURITY™ C₁₈

16



Ideal for LC/MS applications.

- Rigorously quality tested throughout manufacturing process
- First-class efficiency and performance
- Highly reproducible, rugged and robust

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-100-020P	22103-154630	4.6	3	150
HSA-115-515U	22105-154630	4.6	5	150
HSA-115-525R	22105-254630	4.6	5	250

Guard columns

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-100-600U	22103-014001	4.0 and 4.6	3	10	4
HSA-115-600G	22105-014001	4.0 and 4.6	5	10	4

HPLC columns, Thermo Scientific HyPURITY™ ADVANCE

16



- A polar embedded stationary phase for alternate selectivity
- Enhanced aqueous stability
- Exceptional peak shape for acids, bases and chelators

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-050-020L	21003-154630	4.6	3	150
HSA-050-050C	21005-154630	4.6	5	150
HSA-050-060W	21005-254630	4.6	5	250

Guard columns

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-050-610N	21003-014001	4.0 and 4.6	3	10	4
HSA-050-600Q	21005-014001	4.0 and 4.6	5	10	4