

Dispensers, bottle top, Varispenser® and Varispenser® Plus

eppendorf

a

97

Technical Specification - Specific

Volume, mL	Subdivisions, mL	Inaccuracy, %	Imprecision, %
0.50 to 2.50	0.05	±0.6	≤0.1
1.00 to 5.00	0.10	±0.5	≤0.1
2.00 to 10.0	0.20	±0.5	≤0.1
5.00 to 25.0	0.50	±0.5	≤0.1
10.0 to 50.0	1.00	±0.5	≤0.1
20.0 to 100.0	2.00	±0.5	≤0.1

Test conditions

Liquid	Bi-distilled water
Reference temperature	20 to 25°C, constant to ±0.5°C
Number of determinations	10, in accordance with DIN 12650

Varispenser®

Catalogue No	Alt. No	Volume range, mL	O.D., mm
DHT-150-010Y	4961000012	0.50 to 2.50	28, 40, 45 (for adapters)
DHT-150-030C	4961000020	1.00 to 5.00	28, 40, 45 (for adapters)
DHT-150-050T	4961000039	2.0 to 10.0	28, 40, 45 (for adapters)
DHT-150-070N	4961000047	5.0 to 25.0	38, 40, 32/45 (for adapters)
DHT-150-090H	4961000055	10.0 to 50.0	38, 40, 32/45 (for adapters)
DHT-150-110E	4961000063	20.0 to 100.0	38, 40, 32/45 (for adapters)

Varispenser®

Catalogue No	Alt. No	Volume range, mL	O.D., mm
DHT-155-010B	4960000019	0.50 to 2.50	28, 40, 45 (for adapters)
DHT-155-030S	4960000027	1.00 to 5.00	28, 40, 45 (for adapters)
DHT-155-050M	4960000035	2.0 to 10.0	28, 40, 45 (for adapters)
DHT-155-070G	4960000043	5.0 to 25.0	38, 40, 32/45 (for adapters)
DHT-155-090A	4960000051	10.0 to 50.0	38, 40, 32/45 (for adapters)
DHT-155-110U	4960000060	20.0 to 100.0	38, 40, 32/45 (for adapters)



Dispensers, bottle top, 'Ultra'

03

volac

a

Bottle top dispenser to dispense set amounts of liquids.

The standard pack includes dispenser unit, 80mL amber reservoir (10mL divisions), with PTFE lined cap and two adapters, GL45 and 32mm.

The spring loaded plunger allows use as a repeater pipettor from a reservoir, or as an alternative to a 'stepper'. Lower body is autoclavable to 121°C.

Catalogue No	Alt. No	Volume, mL	Subdivisions, µL	Accuracy, %
DHT-256-010K	AR675/O	0.05 to 0.5	1	1
DHT-256-030E	AR675/A	0.1 to 1	2	1
DHT-256-050V	AR675/B	0.5 to 5	10	1
DHT-256-070P	AR675/C	1 to 10	20	1

Catalogue No	Alt. No	Description	Capacity, mL
DHT-258-050J	RSMD16	Reservoir for dispenser	400
DHT-258-070D	RSMD17	Reservoir for dispenser	80

