





Filters, membrane, PTFE, Fluoropore™ and Mitex™ 71



Fluoropore™ - Hydrophobic PTFE membrane bonded to a high density polyethylene support. Clarification of acids, bases and solvents, clarification of alcohol solutions, UV spectroscopy. Wetting with methanol or another low surface tension fluid renders the membrane hydrophilic, mean thickness varies from 145µm to 200µm.

Mitex™ - An unsupported hydrophobic PTFE membrane. Clarification of acids, bases and cryogenic fluids, clarification of propellants, analysis of hydraulic fluids, RNA isolation. Mean thickness is 125µm, withstands temperatures up to 260°C.

Autoclavable at 121°C

Hydrophobic PTFE Millipore membranes and LCR and Omnipore membranes are also available. Please contact Customer Service.

Fluoropore™ membrane filters
0.2µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-330-010A	FGLP01300	13	Plain	100
FDR-330-020U	FGLP 025 00	25	Plain	100
FDR-330-030R	FGLP04700	47	Plain	100
FDR-330-040X	FGLP09050	90	Plain	50
FDR-330-050L	FGLP14250	142	Plain	50

0.5µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-331-010R	FHLP01300	13	Plain	100
FDR-331-020X	FHLP02500	25	Plain	100
FDR-331-030L	FHLP03700	37	Plain	100
FDR-331-040Y	FHLP04700	47	Plain	100
FDR-331-060C	FHLP14250	142	Plain	50

0.5µm pore size, white, unlaminated

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-332-010L	FHUP04700	47	Plain	100

1.0µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-333-010F	FALP01300	13	Plain	100
FDR-333-020C	FALP02500	25	Plain	100
FDR-333-030W	FALP04700	47	Plain	100
FDR-333-040T	FALP09050	90	Plain	50
FDR-333-050Q	FALP14250	142	Plain	50

3.0µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-334-010W	FSLW02500	25	Plain	100
FDR-334-020T	FSLW04700	47	Plain	100
FDR-334-030Q	FSLW09025	90	Plain	25
FDR-334-040N	FSLW14200	142	Plain	10

Mitex™ membrane filters

5.0µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-335-010Q	LSWP01300	13	Plain	100
FDR-335-020N	LSWP02500	25	Plain	100
FDR-335-030K	LSWG02500	25	Grid	100
FDR-335-040H	LSWP03700	37	Plain	100
FDR-335-050E	LSWP04700	47	Plain	100
FDR-335-060B	LSWG04700	47	Grid	100
FDR-335-070V	LSWP09025	90	Plain	25
FDR-335-080S	LSWP14250	142	Plain	50

10.0µm pore size, white

Catalogue No	Alt. No	Diameter, mm	Type	Pack qty
FDR-336-010K	LCWP01300	13	Plain	100
FDR-336-020H	LCWP02500	25	Plain	100
FDR-336-030E	LCWG02500	25	Grid	100
FDR-336-040B	LCWP04700	47	Plain	100
FDR-336-050V	LCWG04700	47	Grid	100
FDR-336-060S	LCWP09025	90	Plain	25
FDR-336-070P	LCWP14250	142	Plain	50

