



## Adapters, multiple neck

01

### 1 Two necks, parallel vertical

Catalogue No	Alt. No	Socket size	Cone size
QAF-020-J	MA1/00	10/19	10/19
QAF-030-C	MA1/11	14/23	14/23
QAF-040-S	MA1/2	19/26	19/26
QAF-050-L	MA1/3	19/26	24/29
QAF-060-E	MA1/33	24/29	24/29

### 2 Two necks, one vertical and one at 45°

Catalogue No	Alt. No	Socket size	Cone size
QAF-200-H	MA2/2	19/26	19/26
QAF-210-A	MA2/3	19/26	24/29

### 3 Three necks, two parallel vertical and one at 45°

Catalogue No	Alt. No	Socket size	Cone size
QAF-410-H	MA4/23	19/26 (parallel) 19/26 (angled)	24/29

## Adapters, receiver

01

### Multi-limb adapter

Allows four, 100mL. flasks to be connected to a condenser outlet. Can be rotated at the conical joint allowing each flask to be progressively filled without disturbing the vacuum within the delivery system.

### Perkin adapters

For vacuum distillation work. Receiver flask can be changed without disturbing the vacuu within the system. Fitted with 3 x GP Rotaflo stopcocks with PTFE keys. With screwcap system for gentle release of vacuum. Without flask.



### Receiver, vacuum, multi-limb

Catalogue No	Alt. No	Socket	Cone size
QAF-501-A	RA11/SC	19/26	14/23

### Vertical with sidearm for vacuum connection

Catalogue No	Alt. No	Socket	Cone size
QAH-591-J	RA13/33SC	24/29	24/29

### Intermediate (Perkin)

Catalogue No	Alt. No	Socket	Cone size	Capacity, mL
QAF-531-C	RA6/22SC	19/26	19/26	32 (minimum)
QAF-541-S	RA6/23SC	19/26	24/29	32 (minimum)
QAF-551-L	RA6/33SC	24/29	24/29	32 (minimum)



## Bellows, flexible, PTFE, for jointed glassware

42

PTFE flexible bellows connectors are machined from pure PTFE and are used to correct misalignment and relieve stress in glassware assemblies. Connect directly to standard ground glass joints where trapezoidal sealing rings give medium to high vacuum performance. Inert and usable to approx 280°C.

Catalogue No	Socket size	Cone size	Height, mm
FB69230	14/23	14/23	90
FB69231	14/23	29/32	110
FB69232	19/26	19/26	100
FB69233	19/26	24/29	105
FB69234	24/29	24/29	110
FB69235	29/32	29/32	115
FB69236	29/32	34/35	120
FB69237	34/35	34/35	120

