



Newton force meter, digital

- Digital force meter with metal hook and built in measuring tape designed to easily measure force in many applications
- Tare function (used to subtract weight of container)
- Automatic shut off
- Capacity: 150N x 0.2N
- Dimensions (w x d x h), mm: 90 x 30 x 270
- Requires 1 x 9V 6LR22 alkaline battery (Not included)

Catalogue No
XBY-450-Q

Accessory
Catalogue No
BMT-680-090H

Description
Alkaline
battery 9V
(6LR22)



Force meter, platform, Newton

- Platform type force meter for personal weighing with zero adjustment and magnified scale window
- Capacity: 1,275N x 5N
- Dimensions (w x d x h), mm: 310 x 245 x 60

Catalogue No
XBY-500-U



Masses, iron, hexagonal, with lifting ring

Catalogue No	Mass
MAS-600-050Q	50g
MAS-600-060N	100g
MAS-600-070K	200g
MAS-600-080H	500g
MAS-600-090E	1kg
MAS-600-100E	2kg
MAS-600-110B	5kg
MAS-600-120V	10kg



Masses, hanger with slotted 100g masses

- Comprises 100g mass hanger with hook, cast alloy base and nine alloy 100g masses, accuracy 5%, total mass 1kg
- Slot prevents masses sliding off hanger

Catalogue No
MAS-700-A

Pack qty
4

Spare Parts

Catalogue No	Description	Pack qty
MAS-704-030N	Masses, 100g	5
MAS-704-050H	Mass hangers, 100g	5



**Balances, General Labware,
pages 365 to 373**



Masses, weighing set on a tray

NEW

- Combined set of accurate steel and brass weights
- Housed in sturdy moulded tray
- For weighing objects from 10g to 1,990g in 10g increments
- 3 steel weights and 6 brass weights
- Contents:
1 each of 1kg, 500g, 200g, 50g, 20g
2 each of 100g, 10g

Catalogue No
MAS-200-010X